

V300 E

ELECTRONIC FOLDABLE WHEELCHAIR

Electronic wheelchair equipped with footrests and armrests adjustable in height and width. Equipped with a powerful engine of 2 x 150W and 26Ah lithium battery.



C76
White

Pushbar

Foldable backrest

Armrests
- Adjustable
- Removable

Legrests
- Adjustable
- Removable



Folded 420(W) x 830(L) x (760-770)(H)



24V-26Ah 2x150W 6 km/h 20 km 10° 15 1400



605	420	420	200 - 300	460	450	970	450	1060	40	90	440(W)
645		460									850(L)
685		500									800(H)

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1.5 kg and 1.5°) are possible.
Specifications might change without prior notice and are therefore to be confirmed.
* Special delivery conditions apply.