

# Express

## ELECTRONIC FOLDABLE WHEELCHAIR

Electronic wheelchair equipped with footrests and armrests adjustable in height and width. Equipped with a powerful engine of 2 x 150W and 38 A/h or 50A/h batteries.



**C59**

Champagne



### Wheelchair Control System

- Ideal for light weight, indoor wheelchairs
- Easy to use
- Reliable



Folded 400(W) x 670(L) x (670-750)(H)



Pushbar

Armrests  
- Adjustable  
- Removable

Legrests  
- Adjustable  
- Removable



2x12V-38Ah	2x150W	6 km/h	18 km	7°	15	1400
*2x12V-50Ah			25 km			



610	400	420	225 · 315	540	430	930	430	950	67	90	650(W)
630		460									850(L)
690		500									650(H)

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1.5 kg and 1.5°) are possible.  
Specifications might change without prior notice and are therefore to be confirmed.  
\* Special delivery conditions apply.